

Human RLIM protein (Recombinant) (N-His) (STJP007928)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	32.94 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	51132
Gene Symbol	RLIM
UniProt ID	RNF12_HUMAN
Immunogen Region	Ser97-Ser370

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--