

Human NDUFB6 protein (Recombinant) (N-His-SUMO & C-Strep) (STJP008119)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli
PRODUCT PROPERTIES	
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	21.86 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	4712
Gene Symbol	NDUFB6
UniProt ID	NDUB6_HUMAN
Immunogen Region	Thr2-Ser68
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.