

EBV/HHV4 gH & gL & gp42 Trimer protein (Recombinant) (C-Strep & C-His & C-His) (STJP008444)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	74 kDa & 14 kDa & 25 kDa (Reduced)
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	gH [Ser19-Leu682] & gL [Asn23-Gly137] & gp42 [Arg30-Ser223]
------------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--