

## Human TUBA1B/Alpha Tubulin protein (Recombinant) (C-His) (STJP009355)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	51.21 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">10376</a>
Gene Symbol	<a href="#">TUBA1B</a>
UniProt ID	<a href="#">TBA1B_HUMAN</a>
Immunogen Region	Met1-Tyr451

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--