

## Human BP230 protein (Recombinant) (STJP009372)

### GENERAL INFORMATION

Product Type    Proteins

### PRODUCT PROPERTIES

Formulation    Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.  
Storage Instruction    Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID    [667](#)  
Gene Symbol    [DST](#)  
UniProt ID    [DYST\\_HUMAN](#)

### ADDITIONAL INFORMATION

Note    STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.