

Mouse KLRG1/CLEC15A protein (Recombinant) (N-GST & C-His) (STJP009778)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	41.12 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	50928
Gene Symbol	Klrp1
UniProt ID	KLRG1_MOUSE
Immunogen Region	Ile77-Tyr188

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--