

**BKPyV Agnoprotein protein (Recombinant) (N-His-SUMO) (STJP010429)****GENERAL INFORMATION**

Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli

**PRODUCT PROPERTIES**

Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	17.69 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

**TARGET INFORMATION**

Immunogen Region	[Met1-Arg27] & [Ser45-Ser66]
------------------	------------------------------

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--