

**Brucella suis omp19 protein (Recombinant) (N-GST & C-His) (STJP010559)**

| GENERAL INFORMATION    |                                                                                                                                                                           |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type           | Proteins                                                                                                                                                                  |
| Applications           | ELISA/Immunogen/SDS-PAGE/WB                                                                                                                                               |
| Host / Source          | E.coli                                                                                                                                                                    |
| PRODUCT PROPERTIES     |                                                                                                                                                                           |
| Formulation            | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                |
| Molecular Weight       | 38.06 kDa                                                                                                                                                                 |
| Storage Instruction    | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |
| TARGET INFORMATION     |                                                                                                                                                                           |
| Immunogen Region       | Pro80-Arg177                                                                                                                                                              |
| ADDITIONAL INFORMATION |                                                                                                                                                                           |
| Note                   | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                                      |