

Human SLC2A10 protein (Recombinant) (N-His-KSI) (STJP010789)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli
PRODUCT PROPERTIES	
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	25.55 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	81031
Gene Symbol	SLC2A10
UniProt ID	GTR10_HUMAN
Immunogen Region	Ala320-Leu414
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.