

Biotin-labeled Human IL6 protein (Recombinant) (C-His) {Biotin} (STJP011499)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	Mammalian cells

PRODUCT PROPERTIES

Formulation	0.01M PBS, pH 7.4, 0.2% BSA, 50% Glycerol, 0.05% Proclin 300.
Molecular Weight	25 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3569
Gene Symbol	IL6
UniProt ID	IL6_HUMAN
Immunogen Region	Val30-Met212

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--