

CA16 P2A/Protease 2A protein (Recombinant) (N-His) (STJP011519)**GENERAL INFORMATION**

Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	18.80 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	Gly863-Gln1012
------------------	----------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--