

## CA16 Protease 3C protein (Recombinant) (N-His) (STJP011523)

| GENERAL INFORMATION    |                                                                                                                                                                       |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type           | Proteins                                                                                                                                                              |
| Applications           | ELISA/Immunogen/SDS-PAGE/WB                                                                                                                                           |
| Host / Source          | E.coli                                                                                                                                                                |
| PRODUCT PROPERTIES     |                                                                                                                                                                       |
| Formulation            | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                            |
| Molecular Weight       | 22.32 kDa                                                                                                                                                             |
| Storage Instruction    | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |
| TARGET INFORMATION     |                                                                                                                                                                       |
| Immunogen Region       | Gly1549-Gln1731                                                                                                                                                       |
| ADDITIONAL INFORMATION |                                                                                                                                                                       |
| Note                   | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                                  |