

Influenza A virus HA/Hemagglutinin protein (Recombinant) (C-8His) (STJP011591)

GENERAL INFORMATION

| | |
|---------------|-----------------------------|
| Product Type | Proteins |
| Applications | ELISA/Immunogen/SDS-PAGE/WB |
| Host / Source | Mammalian cells |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Formulation | Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol. |
| Molecular Weight | 60.52 kDa |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|--------------|
| Immunogen Region | Gln17-Asp529 |
|------------------|--------------|

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|