

Influenza B virus HA/Hemagglutinin protein (Recombinant) (C-6His) (STJP011614)

GENERAL INFORMATION		
	Product Type	Proteins
	Applications	ELISA/Immunogen/SDS-PAGE/WB
	Host / Source	Mammalian cells
PRODUCT PROPERTIES		
	Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
	Molecular Weight	60.30 kDa
	Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION		
	Immunogen Region	Asp16-Ile548
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.