

CU43 Protease (Human IgG Specific) protein (Recombinant) (STJP011631)**GENERAL INFORMATION**

Product Type	Proteins
Host / Source	Mammalian cells

PRODUCT PROPERTIES

Formulation	0.01M PBS, pH 7.4.
Molecular Weight	70.94 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	CU43 Protease is cloned from <i>Corynebacterium ulcerans</i> and expressed in Mammalian cells.
------------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--