

Mouse SEMA3A protein (Recombinant-Active) (C-hFc-His) (STJP011717)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells
PRODUCT PROPERTIES	
Formulation	Lyophilized from a solution in PBS pH 7.4, 5% Trehalose, 5% Mannitol.
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	20346
Gene Symbol	Sema3a
UniProt ID	SEM3A_MOUSE
Immunogen Region	Asn21-Lys747 (Arg552Ala, Arg554Ala, Arg555Ala, Arg618Ala, Lys619Ala, Lys642Ala, Lys643Ala)
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.