

Human CD256/TNFSF13/APRIL protein (Recombinant-Active) (N-His) (STJP011748)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 5% Trehalose, 5% Mannitol.
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	8741
Gene Symbol	TNFSF13
UniProt ID	TNF13_HUMAN
Immunogen Region	Lys112-Leu250

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--