

Human CD336/NCR2 protein (Recombinant-Active) (C-His) (STJP011759)

GENERAL INFORMATION

| | |
|---------------|-----------------|
| Product Type | Proteins |
| Applications | ELISA/WB |
| Host / Source | Mammalian Cells |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Formulation | Lyophilized from a solution in PBS pH 7.4, 5% Trehalose, 5% Mannitol. |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8Å°C for one week. Store at -20 to -80Å°C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene ID | 9436 |
| Gene Symbol | NCR2 |
| UniProt ID | NCTR2_HUMAN |
| Immunogen Region | Gln22-Pro190 |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|