

Human IDS/Idursulfase protein (Recombinant-Active) (C-10*His) (STJP011926)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 5% Trehalose, 5% Mannitol.
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3423
Gene Symbol	IDS
UniProt ID	IDS_HUMAN
Immunogen Region	Ser26-Pro550

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--