

Human LY6E/SCA-2 protein (Recombinant-Active) (C-hFc-His) (STJP012122)

GENERAL INFORMATION

| | |
|---------------|-----------------|
| Product Type | Proteins |
| Applications | ELISA/WB |
| Host / Source | Mammalian Cells |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Formulation | Lyophilized from a solution in PBS pH 7.4, 5% Trehalose, 5% Mannitol. |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|----------------------------|
| Gene ID | 4061 |
| Gene Symbol | LY6E |
| UniProt ID | LY6E_HUMAN |
| Immunogen Region | Leu21-Ser101 |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|