

## Mouse RSP01 protein (Recombinant-Active) (C-His) (STJP012334)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | ELISA/WB |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Formulation         | Lyophilized from a solution in PBS pH 7.4, 5% Trehalose, 5% Mannitol.   |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

### TARGET INFORMATION

|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">192199</a>      |
| Gene Symbol      | <a href="#">Rspo1</a>       |
| UniProt ID       | <a href="#">RSP01_MOUSE</a> |
| Immunogen Region | Ser21-Gly209                |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|