

Human Laminin 511 E8 protein (Recombinant-Active) (N-His & N-Flag & N-HA) (STJP012353)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.2 µm filtered solution in PBS, pH 7.4.
Molecular Weight	86.7 kDa & 26.4 kDa & 28.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	[Ala2534-Ala3323] & [Leu1561-Leu1786] & [Asn1364-Pro1609]
------------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--