

Human ANGPTL7 protein (Recombinant-Active) (C-10*His) (STJP012357)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mg/ml filtered solution in 20mM PB, 6% Trehalose, 4% Mannitol, 100mM NaCl, 0.05% Tween 80, pH 6.0.
Molecular Weight	38.4 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	10218
Gene Symbol	ANGPTL7
UniProt ID	ANGL7_HUMAN
Immunogen Region	Gln27-Pro346

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--