

Human FGF17 protein (Recombinant-Active) (C-His) (STJP012359)

GENERAL INFORMATION

| | |
|---------------|-----------------|
| Product Type | Proteins |
| Applications | ELISA/WB |
| Host / Source | Mammalian Cells |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Formulation | Lyophilized from a 0.2 Mu m filtered solution in PBS, pH 7.4. |
| Molecular Weight | 22.64 kDa |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene ID | 8822 |
| Gene Symbol | FGF17 |
| UniProt ID | FGF17_HUMAN |
| Immunogen Region | Thr23-Thr216 |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|