

Human PDGFB protein (Recombinant-Active) (Tag free) (STJP012377)

GENERAL INFORMATION

| | |
|---------------|----------|
| Product Type | Proteins |
| Applications | ELISA/WB |
| Host / Source | E.coli |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Formulation | Lyophilized from a 0.2 Mu m filtered solution in 20 mM NaAc-HAc, pH 4.5. |
| Molecular Weight | 12.42 kDa |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene ID | 5155 |
| Gene Symbol | PDGFB |
| UniProt ID | PDGFB_HUMAN |
| Immunogen Region | Ser82-Thr190 |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|