

Biotinylated Human TGFB1/TGF-beta-1 protein (Active) (N-Avi) {Biotin} (STJP012381)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.2 Mu m filtered solution in 50 mM Glycine-HCl, 150 mM NaCl, pH 2.5.
Molecular Weight	14.6 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	7040
Gene Symbol	TGFB1
UniProt ID	TGFB1_HUMAN
Immunogen Region	Ala279-Ser390

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--