

Human INS/Insulin (Active) protein (Recombinant-Active) (STJP012386)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Yeast

PRODUCT PROPERTIES

Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
---------------------	---

TARGET INFORMATION

Gene ID	3630
Gene Symbol	INS
UniProt ID	INS_HUMAN
Immunogen Region	NA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--