

Human IFNB1/IFN-beta protein (Recombinant-Active) (N-Fc) (STJP012393)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.2 µm filtered solution in 20 mM NaAC-HAC, 10%Terhalose, 50 mM NaCl, 3%Mannitol, 0.01%Tween20, pH5.5.
Molecular Weight	46.2kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3456
Gene Symbol	IFNB1
UniProt ID	IFNB_HUMAN
Immunogen Region	Met22-Asn187

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--