

Human IL1B/IL1F2 protein (Recombinant-Active) (Tag free) (STJP012396)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.2 Mu m filtered solution in 20 mM PB, 150 mM NaCl, 3% Trehalose, 0.02% Tween 80, 4% Mannitol, pH 7.4.
Molecular Weight	17.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3553
Gene Symbol	IL1B
UniProt ID	IL1B_HUMAN
Immunogen Region	Ala117-Ser269

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	----------------------------------------------------------------------------------------------------------------------