

Human IGF1 protein (Recombinant-Active) (Tag free) (STJP012409)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.2 Mu m filtered solution in 20 mM NaAc-HAc, pH 4.5.
Molecular Weight	7.3 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3479
Gene Symbol	IGF1
UniProt ID	IGF1_HUMAN
Immunogen Region	Thr52-Ala118

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--