

## Human IL4 protein (Recombinant-Active) (Tag free) (STJP012411)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.2 µm filtered solution in 20 mM Tris-HCl, 150 mM NaCl, pH 7.2.
Molecular Weight	15.1 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">3565</a>
Gene Symbol	<a href="#">IL4</a>
UniProt ID	<a href="#">IL4_HUMAN</a>
Immunogen Region	His25-Ser153

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--