

Human FGF4 protein (Recombinant-Active) (Tag free) (STJP012435)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.2 µm filtered solution in 20 mM PB, 5% Sucrose, 4% Mannitol, 0.02% Tween 80, pH 7.4.
Molecular Weight	16.9 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	2249
Gene Symbol	FGF4
UniProt ID	FGF4_HUMAN
Immunogen Region	Ser54-Leu206

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--