

## Human FGF2/FGFb protein (Recombinant-Active) (Tag free) (STJP012444)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Supplied as a 0.2 µg/ml filtered solution in 20 mM Citrate, 10% Trehalose, 150 mM NaCl, 0.04% PS80, 0.5 mM EDTA, pH 5.5.
Molecular Weight	17.2 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Immunogen Region	Met1-Ser155
------------------	-------------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--