

Human Tau-441 (2N4R) protein (Recombinant-Active) (Tag free) (STJP012459)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mg/ml filtered solution in 20mM PB, 150mM NaCl, 1mM EDTA, pH 7.4.
Molecular Weight	45.85 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	Met1-Leu441
------------------	-------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--