

Human BMP4 (Active) protein (Recombinant-Active) (STJP012467)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.2 μ g/ml filtered solution of 4mM HCl.
Molecular Weight	13.3 KDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at $\leq -80^{\circ}\text{C}$, stable for 1 year after receipt. Store at $\leq -80^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.

TARGET INFORMATION

Gene ID	652
Gene Symbol	BMP4
UniProt ID	BMP4_HUMAN
Immunogen Region	Ser293-Arg408

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--