

Human IL7 protein (Recombinant-Active) (Tag free) (STJP012469)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.2 µm filtered solution in 20 mM PB, 300 mM NaCl, pH 8.0.
Molecular Weight	17.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3574
Gene Symbol	IL7
UniProt ID	IL7_HUMAN
Immunogen Region	Asp26-His177

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--