

Human HGF protein (Recombinant-Active) (C-His) (STJP012472)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.2 µm filtered solution in 20 mM Tris, 500 mM NaCl, 3% Trehalose, 0.02% Tween 80, pH 8.0.
Molecular Weight	26 kDa & 53.7 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3082
Gene Symbol	HGF
UniProt ID	HGF_HUMAN
Immunogen Region	Gln32-Ser728

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--