

## Human LIF protein (Recombinant-Active) (Tag free) (STJP012477)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.2 µm filtered solution in 20 mM PB, 2% Sucrose, 4% Mannitol, 0.02% Tween 80, pH 7.4.
Molecular Weight	19.7 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">3976</a>
Gene Symbol	<a href="#">LIF</a>
UniProt ID	<a href="#">LIF_HUMAN</a>
Immunogen Region	Ser23-Phe202

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--