

Human CD152/CTLA4 protein (Recombinant-Active) (C-Flag) (STJP012494)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM PB, 50mM NaCl, 8% Trehalose, 4% Mannitol, 0.02% Tween 80, pH6.6
Molecular Weight	14.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	1493
Gene Symbol	CTLA4
UniProt ID	CTLA4_HUMAN
Immunogen Region	Lys36-Asp161

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--