

Human LGALS3/Galectin-3 protein (Recombinant-Active) (Tag free) (STJP012497)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 50mM HEPES, 5% Sucrose, 5% Mannitol, 0.06% Tween 80, pH 7.5.
Molecular Weight	26 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3958
Gene Symbol	LGALS3
UniProt ID	LEG3_HUMAN
Immunogen Region	Ala2-Ile250

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--