

Human BMP7 protein (Recombinant-Active) (tag free) (STJP012499)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.2 µg/ml filtered solution in 20 mM PB, 100 mM NaCl, 5% Mannitol, 0.05% Tween 80, pH 7.5.
Molecular Weight	15.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	655
Gene Symbol	BMP7
UniProt ID	BMP7_HUMAN
Immunogen Region	Ser293-His431

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--