

## Biotinylated Human IGFBP5 protein (C-His-Avi) {Biotin} (STJP012538)

### GENERAL INFORMATION

|               |                 |
|---------------|-----------------|
| Product Type  | Proteins        |
| Applications  | ELISA/WB        |
| Host / Source | Mammalian Cells |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilized from a 0.22 Mu m filtered solution in 20mM PB, 300mM NaCl, pH 7.4.                                                                                        |
| Molecular Weight    | 31.4 kDa                                                                                                                                                              |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

### TARGET INFORMATION

|                  |                            |
|------------------|----------------------------|
| Gene ID          | <a href="#">3488</a>       |
| Gene Symbol      | <a href="#">IGFBP5</a>     |
| UniProt ID       | <a href="#">IBP5_HUMAN</a> |
| Immunogen Region | Leu21-Glu272               |

### ADDITIONAL INFORMATION

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|