

Human EPHA8 protein (Recombinant-Active) (C-His) (STJP012556)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM Citrate, 6% Trehalose, 50mM NaCl, 0.05% Tween 80, pH 5.0.
Molecular Weight	57.4 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	2046
Gene Symbol	EPHA8
UniProt ID	EPHA8_HUMAN
Immunogen Region	Glu31-Thr542

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--