

Human CD154/CD40LG/TNFSF5 protein (Recombinant-Active) (Tag free) (STJP012561)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mg/ml filtered solution in 20mM PB, 200mM NaCl, 0.1mM EDTA, pH 7.0.
Molecular Weight	16.2 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	959
Gene Symbol	CD40LG
UniProt ID	CD40L_HUMAN
Immunogen Region	Met113-Leu261

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--