

Human TNFSF9/4-1BBL protein (Recombinant-Active) (N-His-Flag) (STJP012589)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 M filtered solution in 20mM Tris-HCl, 50mM NaCl, 5% Trehalose, 5% Mannitol, 1mM EDTA, pH 8.0.
Molecular Weight	24.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	8744
Gene Symbol	TNFSF9
UniProt ID	TNFL9_HUMAN
Immunogen Region	Arg71-Glu254

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--