

Human CD253/TNFSF10/TRAIL protein (Recombinant) (Tag free) (STJP012606)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in 40mM Tris-HCl, 300mM NaCl, 5% Trehalose, 5% Mannitol, 0.01% Tween80, 10% Glycerol, pH 7.0.
Molecular Weight	19.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	8743
Gene Symbol	TNFSF10
UniProt ID	TNF10_HUMAN
Immunogen Region	Arg115-Gly281

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--