

Human INHBC protein (Recombinant-Active) (N-His) (STJP012612)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM Citrate, 50mM NaCl, 8% Trehalose, 2% Mannitol, 0.02% Tween 80, 1mM DTT, pH4.1
Molecular Weight	14.83 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3626
Gene Symbol	INHBC
UniProt ID	INHBC_HUMAN
Immunogen Region	Gly237-Ser352

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--