

## Human IL2 (aldesleukin) protein (Recombinant-Active) (Tag free) (STJP012619)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | ELISA/WB |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Formulation         | Lyophilized from a 0.2 µm filtered solution in 10 mM Acetata-Na, 5% Trehaiose, pH 4.5.  |
| Molecular Weight    | 15.5 kDa  |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

### TARGET INFORMATION

|                  |                           |
|------------------|---------------------------|
| Gene ID          | <a href="#">3558</a>      |
| Gene Symbol      | <a href="#">IL2</a>       |
| UniProt ID       | <a href="#">IL2_HUMAN</a> |
| Immunogen Region | Pro22-Thr153 (Cys145Ser)  |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|