

Human IL2 Superkine (Mutant) protein (Recombinant-Active) (C-His) (STJP012620)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.2 Mu m filtered solution in 4 mM HCl.
Molecular Weight	16.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3558
Gene Symbol	IL2
UniProt ID	IL2_HUMAN
Immunogen Region	Ala21-Thr153 (L100F, R101D, L105V, I106V, I112F)

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--