

Biotinylated Human CD172a/SIRPA protein (C-His-Avi) {Biotin} (STJP012626)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 µg/ml filtered solution in 20mM PB, 150mM NaCl, pH 7.4.
Molecular Weight	39.9 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	140885
Gene Symbol	SIRPA
UniProt ID	SHPS1_HUMAN
Immunogen Region	Glu31-Arg370

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--